

For Research Use Only

PICALM Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP02979-2

Capture Antibody Information

Catalog Number:
87344-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
252594A3
Reactivity:
human
GenBank:
NM_007166.3
Immunogen Catalog Number:
Eg5966

Conjugate:
Unconjugated
Full name:
phosphatidylinositol binding clathrin assembly protein
Gene ID:
8301

Detection Antibody Information

Catalog Number:
87344-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
252594B7
Reactivity:
human
GenBank:
NM_007166.3
Immunogen Catalog Number:
Eg5966

Conjugate:
Unconjugated
Full name:
phosphatidylinositol binding clathrin assembly protein
Gene ID:
8301

Applications

Tested Applications:
Sandwich ELISA

Range:
0.195-12.5 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02979-2 targets PICALM in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: PICAL Recombinant monoclonal antibody, PBS Only (Capture) 87344-5-PBS (252594A3). 100 µg. Concentration 1 mg/mL.

Detection antibody: PICAL Recombinant monoclonal antibody, PBS Only (Detector) 87344-3-PBS (252594B7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02979-2, Human PICALM Recombinant Matched Antibody Pair - PBS only. 87344-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg5966. 87344-3-PBS was HRP conjugated as the detection antibody. Range: 0.195-12.5 ng/mL



The mean PICALM concentration was determined to be 63.00 ng/mL in K-562 cell extract based on a 1.20 mg/mL extract load, 31.01 ng/mL in HeLa cell extract based on a 1.20 mg/mL extract load and 26.77 ng/mL in HCT 116 cell extract based on a 1.30 mg/mL extract load.