

For Research Use Only

# ST6GAL1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP02980-1

## Capture Antibody Information

**Catalog Number:**  
87412-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
251718E4  
**Reactivity:**  
human  
**GenBank:**  
BC040009  
**Immunogen Catalog Number:**  
Eg5012

**Conjugate:**  
Unconjugated  
**Full name:**  
ST6 beta-galactosamide alpha-2,6-sialyltransferase 1  
**Gene ID:**  
6480

## Detection Antibody Information

**Catalog Number:**  
87412-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
251718B9  
**Reactivity:**  
human  
**GenBank:**  
BC040009  
**Immunogen Catalog Number:**  
Eg5012

**Conjugate:**  
Unconjugated  
**Full name:**  
ST6 beta-galactosamide alpha-2,6-sialyltransferase 1  
**Gene ID:**  
6480

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-80 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP02980-1 targets ST6GAL1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ST6GAL1 Recombinant monoclonal antibody, PBS Only (Capture) 87412-2-PBS (251718E4). 100 µg. Concentration 1 mg/mL.

Detection antibody: ST6GAL1 Recombinant monoclonal antibody, PBS Only (Detector) 87412-1-PBS (251718B9). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MPO2980-1, ST6GAL1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87412-2-PBS. Detection antibody: 87412-1-PBS. Standard: Eg5012. Range: 0.625-80 ng/mL



The mean ST6GAL1 concentration was determined to be 4,121.9 pg/mL in Human plasma.