

GPR183 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP03015-1**

Capture Antibody Information

Catalog Number:

87492-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

252701A7

Reactivity:

human

GenBank:

BC020752

Conjugate:

Unconjugated

Full name:

G protein-coupled receptor 183

Gene ID:

1880

Detection Antibody Information

Catalog Number:

87492-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

252701B2

Reactivity:

human

GenBank:

BC020752

Conjugate:

Unconjugated

Full name:

G protein-coupled receptor 183

Gene ID:

1880

Applications

Tested Applications:

Cytometric bead array

Range:

1.562-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03015-1 targets GPR183 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GPR183 Recombinant monoclonal antibody, PBS Only (Capture) 87492-2-PBS (252701A7). 100 µg. Concentration 1 mg/mL.

Detection antibody: GPR183 Recombinant monoclonal antibody, PBS Only (Detector) 87492-1-PBS (252701B2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C


Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP03015-1, GPR183 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87492-2-PBS. Detection antibody: 87492-1-PBS. Standard: Ag3035. Range: 1.562-200 ng/mL.



The mean GPR183 concentration was determined to be 77.7 ng/mL in HL-60 cell extract based on a 1.2 mg/mL extract load, 42.5 ng/mL in Raji cell extract based on a 1.2 mg/mL extract load, 35.1 ng/mL in SH-SY5Y cell extract based on a 1.2 mg/mL extract load.