

For Research Use Only

Mouse CD11b Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP03033-2

Capture Antibody Information

Catalog Number:
87419-6-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
252750C3
Reactivity:
mouse
GenBank:
Immunogen Catalog Number:
Eg10463

Conjugate:
Unconjugated
Full name:
integrin alpha M
Gene ID:
16409

Detection Antibody Information

Catalog Number:
87419-7-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
252750E6
Reactivity:
mouse
GenBank:
Immunogen Catalog Number:
Eg10463

Conjugate:
Unconjugated
Full name:
integrin alpha M
Gene ID:
16409

Applications

Tested Applications:
Sandwich ELISA

Range:
39-2500 pg/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03033-2 targets CD11b in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: Mouse ITAM Recombinant monoclonal antibody, PBS Only (Capture) 87419-6-PBS (252750C3). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse ITAM Recombinant monoclonal antibody, PBS Only (Detector) 87419-7-PBS (252750E6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

