

For Research Use Only

Mouse BCAM Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP03036-1

Capture Antibody Information

Catalog Number:
87430-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
251702G10

Reactivity:
mouse

GenBank:
NM_020486.2

Immunogen Catalog Number:
Eg6306

Conjugate:
Unconjugated

Full name:
basal cell adhesion molecule

Gene ID:
57278

Detection Antibody Information

Catalog Number:
87430-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
251702C4

Reactivity:
mouse

GenBank:
NM_020486.2

Immunogen Catalog Number:
Eg6306

Conjugate:
Unconjugated

Full name:
basal cell adhesion molecule

Gene ID:
57278

Applications

Tested Applications:
Cytometric bead array

Range:
0.312-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03036-1 targets BCAM in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse CD239/BCAM Recombinant monoclonal antibody, PBS Only (Capture) 87430-3-PBS (251702G10). 100 μ g. Concentration 1 mg/ml.

Detection antibody: Mouse CD239/BCAM Recombinant monoclonal antibody, PBS Only (Detector) 87430-2-PBS (251702C4). 100 μ g. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP03036-1, MOUSE BCAM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87430-3-PBS. Detection antibody: 87430-2-PBS. Standard: Eg6306. Range: 0.312-20 ng/mL



The mean BCAM concentration was determined to be 1420.5 ng/mL in Mouse kidney tissue extract based on a 7.8 mg/mL extract load, 192.3 ng/mL in Mouse skeletal muscle tissue extract based on a 4.1 mg/mL extract load, 106.5 ng/mL in Mouse brain tissue extract based on a 3.2 mg/mL extract load.