

For Research Use Only



ABCC6 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP03056-2

Capture Antibody Information

Catalog Number:
87456-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
252956G12
Reactivity:
human
GenBank:
BC131732

Conjugate:
Unconjugated
Full name:
ATP-binding cassette, sub-family C (CFTR/MRP), member 6
Gene ID:
368

Detection Antibody Information

Catalog Number:
87456-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
252956H6
Reactivity:
human
GenBank:
BC131732

Conjugate:
Unconjugated
Full name:
ATP-binding cassette, sub-family C (CFTR/MRP), member 6
Gene ID:
368

Applications

Tested Applications:
Sandwich ELISA

Range:
0.625-10 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03056-2 targets ABCC6 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: ABCC6 Recombinant monoclonal antibody, PBS Only (Capture) 87456-5-PBS (252956G12). 100 µg. Concentration 1 mg/mL.

Detection antibody: ABCC6 Recombinant monoclonal antibody, PBS Only (Detector) 87456-4-PBS (252956H6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP03056-2, Human ABCC6 Recombinant Matched Antibody Pair - PBS only. 87456-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard SY03023. 87456-4-PBS was HRP conjugated as the detection antibody. Range: 0.625-10 ng/mL



The mean ABCC6 concentration was determined to be 92.00 ng/mL in mouse kidney tissue extract based on a 8.4 mg/mL extract load and 12.50 ng/mL in mouse brain tissue extract based on a 6.4 mg/mL extract load.