

For Research Use Only

SNAP23 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP03062-1

Capture Antibody Information

Catalog Number:
87283-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
252538C11
Reactivity:
human
GenBank:
NM_003825.4
Immunogen Catalog Number:
Eg6079

Conjugate:
Unconjugated
Full name:
synaptosomal-associated protein,
23kDa
Gene ID:
8773

Detection Antibody Information

Catalog Number:
87283-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
252538A4
Reactivity:
human
GenBank:
NM_003825.4
Immunogen Catalog Number:
Eg6079

Conjugate:
Unconjugated
Full name:
synaptosomal-associated protein,
23kDa
Gene ID:
8773

Applications

Tested Applications:
Sandwich ELISA

Range:
78.1-5000 pg/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03062-1 targets SNAP23 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: SNAP23 Recombinant monoclonal antibody, PBS Only (Capture) 87283-2-PBS (252538C11). 100 μ g. Concentration 1 mg/mL.

Detection antibody: SNAP23 Recombinant monoclonal antibody, PBS Only (Detector) 87283-3-PBS (252538A4). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptgcn.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP03062-1, Human SNAP23 Recombinant Matched Antibody Pair - PBS only. 87283-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg6079. 87283-3-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL

The mean SNAP23 concentration was determined to be 220.5 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load and 104.5 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load.