

For Research Use Only

Mouse CD3 Delta Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP03077-1

Capture Antibody Information

Catalog Number:
87532-8-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
252849H6
Reactivity:
mouse
GenBank:
NM_013487.3
Immunogen Catalog Number:
Eg4665

Conjugate:
Unconjugated
Full name:
CD3 antigen, delta polypeptide
Gene ID:
12500

Detection Antibody Information

Catalog Number:
87532-7-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
252849A1
Reactivity:
mouse
GenBank:
NM_013487.3
Immunogen Catalog Number:
Eg4665

Conjugate:
Unconjugated
Full name:
CD3 antigen, delta polypeptide
Gene ID:
12500

Applications

Tested Applications:
Cytometric bead array

Range:
0.312-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03077-1 targets CD3 Delta in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse CD3D Recombinant monoclonal antibody, PBS Only (Capture) 87532-8-PBS (252849H6). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse CD3D Recombinant monoclonal antibody, PBS Only (Detector) 87532-7-PBS (252849A1). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.



Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptgcn.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP03077-1, MOUSE CD3 Delta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87532-8-PBS. Detection antibody: 87532-7-PBS. Standard: Eg4665. Range: 0.312-40 ng/mL

The mean CD3 Delta concentration was determined to be 16.9 ng/mL in CTLL-2 cell extract based on a 1.2 mg/mL extract load, 4.7 ng/mL in mouse thymus tissue extract based on a 1.2 mg/mL extract load, mouse lung tissue extract were served as a negative control.