

For Research Use Only

RIC8B Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP03094-1

Capture Antibody Information

Catalog Number:
87524-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
252908B11

Reactivity:
human

GenBank:
BC034220

Immunogen Catalog Number:
Ag11888

Conjugate:
Unconjugated

Full name:
resistance to inhibitors of
cholinesterase 8 homolog B (C.
elegans)

Gene ID:
55188

Detection Antibody Information

Catalog Number:
87524-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
252908A7

Reactivity:
human, mouse, rat, pig

GenBank:
BC034220

Immunogen Catalog Number:
Ag11888

Conjugate:
Unconjugated

Full name:
resistance to inhibitors of
cholinesterase 8 homolog B (C.
elegans)

Gene ID:
55188

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP03094-1 targets RIC8B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RIC8B Recombinant monoclonal antibody, PBS Only (Capture) 87524-2-PBS (252908B11). 100 μ g. Concentration 1 mg/mL.

Detection antibody: RIC8B Recombinant monoclonal antibody, PBS Only (Detector) 87524-1-PBS (252908A7). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP03094-1, RIC8B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87524-2-PBS. Detection antibody: 87524-1-PBS. Standard: Ag11888. Range: 0.156-20 ng/mL

The mean RIC8B concentration was determined to be 1.4 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 3.9 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load.