

For Research Use Only

AARS Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP03098-1

Capture Antibody Information

Catalog Number:
87581-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
253004F12

Reactivity:
human

GenBank:
BC011451

Immunogen Catalog Number:
Ag11294

Conjugate:
Unconjugated

Full name:
alanyl-tRNA synthetase

Gene ID:
16

Detection Antibody Information

Catalog Number:
87581-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
253004E4

Reactivity:
human

GenBank:
BC011451

Immunogen Catalog Number:
Ag11294

Conjugate:
Unconjugated

Full name:
alanyl-tRNA synthetase

Gene ID:
16

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03098-1 targets AARS in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: AARS Recombinant monoclonal antibody, PBS Only (Capture) 87581-4-PBS (253004F12). 100 µg. Concentration 1 mg/mL.

Detection antibody: AARS Recombinant monoclonal antibody, PBS Only (Detector) 87581-3-PBS (253004E4). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP03098-1, AARS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87581-4-PBS. Detection antibody: 87581-3-PBS. Standard: Ag11294. Range: 0.156-20 ng/mL

The mean AARS concentration was determined to be 65.3 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 30.7 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 48.3 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load.