

For Research Use Only

Mouse Podocalyxin Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP03126-1

Capture Antibody Information

Catalog Number:
87582-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
252967D3
Reactivity:
mouse
GenBank:
NM_013723.3
Immunogen Catalog Number:
Eg7526

Conjugate:
Unconjugated
Full name:
podocalyxin-like
Gene ID:
27205

Detection Antibody Information

Catalog Number:
87582-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
252967A3
Reactivity:
mouse
GenBank:
NM_013723.3
Immunogen Catalog Number:
Eg7526

Conjugate:
Unconjugated
Full name:
podocalyxin-like
Gene ID:
27205

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03126-1 targets Podocalyxin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Podocalyxin Recombinant monoclonal antibody, PBS Only (Capture) 87582-3-PBS (252967D3). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse Podocalyxin Recombinant monoclonal antibody, PBS Only (Detector) 87582-2-PBS (252967A3). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptgcn.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP03126-1, MOUSE Podocalyxin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87582-3-PBS. Detection antibody: 87582-2-PBS. Standard: Eg7526. Range: 0.156-20 ng/mL

The mean Podocalyxin concentration was determined to be 60.2 ng/mL in mouse kidney tissue extract based on a 8.4 mg/mL extract load, 68.3 ng/mL in mouse lung tissue extract based on a 5.0 mg/mL extract load, 3.2 ng/mL in mouse skeletal muscles tissue extract based on a 4.2 mg/mL extract load.