

MEF2A Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP03216-1**

Capture Antibody Information

Catalog Number: 87660-3-PBS	Clone ID: 253110E4	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: myocyte enhancer factor 2A
Isotype: IgG	GenBank: BC013437	Gene ID: 4205
Purification Method: Protein A purification	Immunogen Catalog Number: Ag41519	

Detection Antibody Information

Catalog Number: 87660-2-PBS	Clone ID: 253110B3	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human, mouse, rat	Full name: myocyte enhancer factor 2A
Isotype: IgG	GenBank: BC013437	Gene ID: 4205
Purification Method: Protein A purification	Immunogen Catalog Number: Ag41519	

Applications

Tested Applications: Sandwich ELISA	Range: 0.156-10 ng/mL (Sandwich ELISA)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	--

Product Information

MP03216-1 targets MEF2A in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: MEF2A Recombinant monoclonal antibody, PBS Only (Capture) 87660-3-PBS (253110E4). 100 µg. Concentration 1 mg/mL.

Detection antibody: MEF2A Recombinant monoclonal antibody, PBS Only (Detector) 87660-2-PBS (253110B3). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only


Selected Validation Data



Sandwich ELISA standard curve of MP03216-1, Human MEF2A Recombinant Matched Antibody Pair - PBS only. 87660-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag41519. 87660-2-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL



Plasma of six individual healthy human donors was measured. The MEF2A concentration of detected samples was determined to be 0.30 ng/mL with a range of 0.07-0.78 ng/mL.



The mean MEF2A concentration was determined to be 4.99 ng/mL in HeLa cell extract based on a 1.20 mg/mL extract load.