

ST3GAL2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP03241-1**

Capture Antibody Information

Catalog Number: 85634-3-PBS	Clone ID: 253479E11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: ST3 beta-galactoside alpha-2,3-sialyltransferase 2
Isotype: IgG	GenBank: BC036777	Gene ID: 6483
Purification Method: Protein A purification	Immunogen Catalog Number: Ag4528	

Detection Antibody Information

Catalog Number: 85634-2-PBS	Clone ID: 253479G5	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: ST3 beta-galactoside alpha-2,3-sialyltransferase 2
Isotype: IgG	GenBank: BC036777	Gene ID: 6483
Purification Method: Protein A purification	Immunogen Catalog Number: Ag4528	

Applications

Tested Applications: Cytometric bead array	Range: 1.562-200 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP03241-1 targets ST3GAL2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ST3GAL2 Recombinant monoclonal antibody, PBS Only (Capture) 85634-3-PBS (253479E11). 100 µg. Concentration 1 mg/mL.

Detection antibody: ST3GAL2 Recombinant monoclonal antibody, PBS Only (Detector) 85634-2-PBS (253479G5). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP03241-1, ST3GAL2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85634-3-PBS. Detection antibody: 85634-2-PBS. Standard: Ag4528. Range: 1.562-200 ng/mL.



The mean ST3GAL2 concentration was determined to be 15.8 ng/mL in HT-29 cell extract based on a 1.2 mg/mL extract load, 0.5 ng/mL in A431 cell extract based on a 1.2 mg/mL extract load.