

SUB1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP03355-2**

Capture Antibody Information

Catalog Number: 88052-3-PBS	Clone ID: 260206G10	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: SUB1 homolog (S. cerevisiae)
Isotype: IgG	GenBank: BC009610	Gene ID: 10923
Purification Method: Protein A purification	Immunogen Catalog Number: Ag1390	

Detection Antibody Information

Catalog Number: 88052-1-PBS	Clone ID: 260206C11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: SUB1 homolog (S. cerevisiae)
Isotype: IgG	GenBank: BC009610	Gene ID: 10923
Purification Method: Protein A purification	Immunogen Catalog Number: Ag1390	

Applications

Tested Applications: Sandwich ELISA	Range: 0.625-20 ng/mL (Sandwich ELISA)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	--	---

Product Information

MP03355-2 targets SUB1 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: SUB1 Recombinant monoclonal antibody, PBS Only (Capture) 88052-3-PBS (260206G10). 100 µg. Concentration 1 mg/mL.

Detection antibody: SUB1 Recombinant monoclonal antibody, PBS Only (Detector) 88052-1-PBS (260206C11). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.



Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Sandwich ELISA standard curve of MP03355-2, Human SUB1 Recombinant Matched Antibody Pair - PBS only. 88052-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag1390. 88052-1-PBS was HRP conjugated as the detection antibody. Range: 0.625-20 ng/mL.

The mean SUB1 concentration was determined to be 523.40 ng/mL in HL-60 cell extract based on a 1.20 mg/mL extract load and 122.10 ng/mL in Ramos cell extract based on a 1.50 mg/mL extract load.