

Syntabulin Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP03456-1**

Capture Antibody Information

Catalog Number: 88360-2-PBS	Clone ID: 260730A9	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: Golgi-localized protein
Isotype: IgG	GenBank: BC012082	Gene ID: 55638
Purification Method: Protein A purification	Immunogen Catalog Number: Ag9344	

Detection Antibody Information

Catalog Number: 88360-1-PBS	Clone ID: 260730D4	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: Golgi-localized protein
Isotype: IgG	GenBank: BC012082	Gene ID: 55638
Purification Method: Protein A purification	Immunogen Catalog Number: Ag9344	

Applications

Tested Applications: Cytometric bead array	Range: 0.312-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP03456-1 targets Syntabulin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GOLSYN Recombinant monoclonal antibody, PBS Only (Capture) 88360-2-PBS (260730A9). 100 µg. Concentration 1 mg/mL.

Detection antibody: GOLSYN Recombinant monoclonal antibody, PBS Only (Detector) 88360-1-PBS (260730D4). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP03456-1, Syntabulin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88360-2-PBS. Detection antibody: 88360-1-PBS. Standard: Ag9344. Range: 0.312-40 ng/mL



The mean Syntabulin concentration was determined to be 4.8 ng/mL in HT-29 cell extract based on a 1.2 mg/mL extract load, 2.5 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load, HL-60 cell extract were served as a negative control.