

For Research Use Only

S100P Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP03518-1

Capture Antibody Information

Catalog Number:
88256-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
260786G1

Reactivity:
human

GenBank:
BC006819

Immunogen Catalog Number:
Ag2379

Conjugate:
Unconjugated

Full name:
S100 calcium binding protein P

Gene ID:
6286

Detection Antibody Information

Catalog Number:
88256-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
260786D7

Reactivity:
human

GenBank:
BC006819

Immunogen Catalog Number:
Ag2379

Conjugate:
Unconjugated

Full name:
S100 calcium binding protein P

Gene ID:
6286

Applications

Tested Applications:
Cytometric bead array

Range:
0.312-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03518-1 targets S100P in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: S100P Recombinant monoclonal antibody, PBS Only (Capture) 88256-1-PBS (260786G1). 100 μ g. Concentration 1 mg/mL.

Detection antibody: S100P Recombinant monoclonal antibody, PBS Only (Detector) 88256-3-PBS (260786D7). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP03518-1, S100P Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88256-1-PBS. Detection antibody: 88256-3-PBS. Standard: Ag2379. Range: 0.312-40 ng/mL



The mean S100P concentration was determined to be 464.9 ng/mL in BxPC-3 cell extract based on a 1.2 mg/mL extract load, SH-SY5Y cell extract were served as a negative control.