

UNC93B1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP03543-2**

Capture Antibody Information

Catalog Number: 88417-1-PBS	Clone ID: 260829D11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: unc-93 homolog B1 (C. elegans)
Isotype: IgG	GenBank: BC101568	Gene ID: 81622
Purification Method: Protein A purification	Immunogen Catalog Number: Ag27937	

Detection Antibody Information

Catalog Number: 88417-3-PBS	Clone ID: 260829A11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: unc-93 homolog B1 (C. elegans)
Isotype: IgG	GenBank: BC101568	Gene ID: 81622
Purification Method: Protein A purification	Immunogen Catalog Number: Ag27937	

Applications

Tested Applications: Sandwich ELISA	Range: 0.391-25 ng/mL (Sandwich ELISA)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	--	---

Product Information

MP03543-2 targets UNC93B1 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: UNC93B1 Recombinant monoclonal antibody, PBS Only (Capture/Detector) 88417-1-PBS (260829D11). 100 µg. Concentration 1 mg/mL.

Detection antibody: UNC93B1 Recombinant monoclonal antibody, PBS Only (Detector) 88417-3-PBS (260829A11). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.



Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Sandwich ELISA standard curve of MP03543-2, Human UNC93B1 Recombinant Matched Antibody Pair - PBS only. 88417-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag27937. 88417-3-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL

The mean UNC93B1 concentration was determined to be 72.87 ng/mL in HL-60 cell extract based on a 1.20 mg/mL extract load and 8.59 ng/mL in IMR-32 cell extract based on a 1.70 mg/mL extract load.