

EIF4G1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP03602-2**

Capture Antibody Information

Catalog Number:
88573-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
261246F9
Reactivity:
human
GenBank:
BC007788
Immunogen Catalog Number:
Ag8342

Conjugate:
Unconjugated
Full name:
eukaryotic translation initiation factor 4 gamma, 1
Gene ID:
1981

Detection Antibody Information

Catalog Number:
88573-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
261246F2
Reactivity:
human
GenBank:
BC007788
Immunogen Catalog Number:
Ag8342

Conjugate:
Unconjugated
Full name:
eukaryotic translation initiation factor 4 gamma, 1
Gene ID:
1981

Applications

Tested Applications:
Sandwich ELISA

Range:
0.391-25 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03602-2 targets EIF4G1 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: EIF4G1 Recombinant monoclonal antibody, PBS Only (Capture) 88573-5-PBS (261246F9). 100 µg. Concentration 1 mg/mL.

Detection antibody: EIF4G1 Recombinant monoclonal antibody, PBS Only (Detector) 88573-4-PBS (261246F2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Sandwich ELISA standard curve of MP03602-2, Human EIF4G1 Recombinant Matched Antibody Pair - PBS only. 88573-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag8342. 88573-4-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL.



The mean EIF4G1 concentration was determined to be 163.05 ng/mL in A431 cell extract based on a 1.20 mg/mL extract load, 127.26 ng/mL in U-251 cell extract based on a 1.20 mg/mL extract load and 123.04 ng/mL in K-562 cell extract based on a 1.20 mg/mL extract load.