

BRAF Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP03686-1

Capture Antibody Information

Catalog Number:
88514-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
260358C10
Reactivity:
human
GenBank:
BC101757
Immunogen Catalog Number:
Ag42145

Conjugate:
Unconjugated
Full name:
v-raf murine sarcoma viral oncogene homolog B1
Gene ID:
673

Detection Antibody Information

Catalog Number:
88514-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
260358A1
Reactivity:
human
GenBank:
BC101757
Immunogen Catalog Number:
Ag42145

Conjugate:
Unconjugated
Full name:
v-raf murine sarcoma viral oncogene homolog B1
Gene ID:
673

Applications

Tested Applications:
Sandwich ELISA

Range:
39.1-2500 pg/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03686-1 targets BRAF in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: BRAF Recombinant monoclonal antibody, PBS Only (Capture) 88514-5-PBS (260358C10). 100 µg. Concentration 1 mg/mL.

Detection antibody: BRAF Recombinant monoclonal antibody, PBS Only (Detector) 88514-4-PBS (260358A1). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Sandwich ELISA standard curve of MP03686-1, Human BRAF Recombinant Matched Antibody Pair - PBS only. 88514-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag42145. 88514-4-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL.



The mean BRAF concentration was determined to be 4,369.8 pg/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 3,460.1 pg/mL in K-562 cell extract based on a 1.2 mg/mL extract load and 1,143.7 pg/mL in U2OS cell extract based on a 1.2 mg/mL extract load.