

NDUFS5 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP03692-1

Capture Antibody Information

Catalog Number:
88827-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
261679H2

Reactivity:
human

GenBank:
BC001884

Immunogen Catalog Number:
Ag7441

Conjugate:
Unconjugated

Full name:
NADH dehydrogenase (ubiquinone) Fe-S protein 5, 15kDa (NADH-coenzyme Q reductase)

Gene ID:
4725

Detection Antibody Information

Catalog Number:
88827-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
261679F6

Reactivity:
human

GenBank:
BC001884

Immunogen Catalog Number:
Ag7441

Conjugate:
Unconjugated

Full name:
NADH dehydrogenase (ubiquinone) Fe-S protein 5, 15kDa (NADH-coenzyme Q reductase)

Gene ID:
4725

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03692-1 targets NDUFS5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NDUFS5 Recombinant monoclonal antibody, PBS Only (Capture) 88827-3-PBS (261679H2). 100 µg. Concentration 1 mg/mL.

Detection antibody: NDUFS5 Recombinant monoclonal antibody, PBS Only (Detector) 88827-2-PBS (261679F6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.


The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP03692-1, NDUFS5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88827-3-PBS. Detection antibody: 88827-2-PBS. Standard: Ag7441. Range: 0.781-100 ng/mL



The mean NDUFS5 concentration was determined to be 225 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 180.3 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load, 126.5 ng/mL in SiHa cell extract based on a 1.2 mg/mL extract load.