

NFATC3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP03718-1**

Capture Antibody Information

Catalog Number:
88911-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
261643H8
Reactivity:
human
GenBank:
BC001050
Immunogen Catalog Number:
Ag12985

Conjugate:
Unconjugated
Full name:
nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 3
Gene ID:
4775

Detection Antibody Information

Catalog Number:
88911-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
261643B6
Reactivity:
human
GenBank:
BC001050
Immunogen Catalog Number:
Ag12985

Conjugate:
Unconjugated
Full name:
nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 3
Gene ID:
4775

Applications

Tested Applications:
Cytometric bead array

Range:
0.078-10 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03718-1 targets NFATC3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NFATC3 Recombinant monoclonal antibody, PBS Only (Capture) 88911-3-PBS (261643H8). 100 µg. Concentration 1 mg/mL.

Detection antibody: NFATC3 Recombinant monoclonal antibody, PBS Only (Detector) 88911-2-PBS (261643B6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP03718-1, NFATC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88911-3-PBS. Detection antibody: 88911-2-PBS. Standard: Ag12985. Range: 0.078-10 ng/mL



The mean NFATC3 concentration was determined to be 4.4 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load, 1.4 ng/mL in PC-3 cell extract based on a 1.2 mg/mL extract load.