

FOXP3 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50216-2

Capture Antibody Information

Catalog Number:
68838-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
3B3E12
Reactivity:
human
GenBank:
BC113401
Immunogen Catalog Number:
Ag17571

Conjugate:
Unconjugated
Full name:
forkhead box P3
Gene ID:
50943

Detection Antibody Information

Catalog Number:
68838-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
4A7E7
Reactivity:
human
GenBank:
BC113401
Immunogen Catalog Number:
Ag17571

Conjugate:
Unconjugated
Full name:
forkhead box P3
Gene ID:
50943

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50216-2 targets FOXP3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FOXP3 Monoclonal antibody, PBS Only (Capture) 68838-2-PBS (3B3E12). 100 µg. Concentration 1 mg/mL.

Detection antibody: FOXP3 Monoclonal antibody, PBS Only (Detector) 68838-3-PBS (4A7E7). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data

