

IL-4RA/CD124 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: MP50242-1

Capture Antibody Information

Catalog Number: 67051-3-PBS	Clone ID: 3F12E8	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: interleukin 4 receptor
Isotype: IgM	GenBank: NM_000418	Gene ID: 3566
Purification Method: Euglobulin precipitation	Immunogen Catalog Number: Ag28411	

Detection Antibody Information

Catalog Number: 67051-2-PBS	Clone ID: 1E9B5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: interleukin 4 receptor
Isotype: IgG1	GenBank: NM_000418	Gene ID: 3566
Purification Method: Protein G purification	Immunogen Catalog Number: Ag28411	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP50242-1 targets IL-4RA/CD124 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IL-4RA/CD124 Monoclonal antibody, PBS Only (Capture/Detector) 67051-3-PBS (3F12E8). 100 µg. Concentration 1 mg/mL.

Detection antibody: IL-4RA/CD124 Monoclonal antibody, PBS Only (Detector) 67051-2-PBS (1E9B5). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50242-1, IL4R/IL-4RA/CD124 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67051-3-PBS. Detection antibody: 67051-2-PBS. Standard:Ag28411. Range: 1.563-100 ng/mL.