

CLCC1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50244-4**

Capture Antibody Information

Catalog Number: 68853-3-PBS	Clone ID: 1F9G8	Conjugate: Unconjugated
Host: Mouse	Reactivity: Human	Full name: chloride channel CLIC-like 1
Isotype: IgG1	GenBank: BC002939	Gene ID: 23155
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag24682	

Detection Antibody Information

Catalog Number: 68853-6-PBS	Clone ID: 3F10C9	Conjugate: Unconjugated
Host: Mouse	Reactivity: Human	Full name: chloride channel CLIC-like 1
Isotype: IgG2b	GenBank: BC002939	Gene ID: 23155
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag24682	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP50244-4 targets CLCC1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CLCC1 Monoclonal antibody, PBS Only (Capture) 68853-3-PBS (1F9G8). 100 µg. Concentration 1 mg/mL.

Detection antibody: CLCC1 Monoclonal antibody, PBS Only (Detector) 68853-6-PBS (3F10C9). 100 µg. Concentration 1 mg/mL.

Alternative CLCC1 matched antibody pairs: MP50244-1, MP50244-2, MP50244-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50244-4, CLCC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68853-3-PBS. Detection antibody: 68853-6-PBS. Standard:Ag24682. Range: 0.781-100 ng/mL.