

STK32A Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50267-1**

Capture Antibody Information

Catalog Number: 67528-1-PBS	Clone ID: 1H2A6	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: serine/threonine kinase 32A
Isotype: IgG2a	GenBank: BC021666	Gene ID: 202374
Purification Method: Protein A purification	Immunogen Catalog Number: Ag10942	

Detection Antibody Information

Catalog Number: 67528-2-PBS	Clone ID: 1B8C7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: serine/threonine kinase 32A
Isotype: IgG1	GenBank: BC021666	Gene ID: 202374
Purification Method: Protein G purification	Immunogen Catalog Number: Ag10942	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-200 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP50267-1 targets STK32A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: STK32A Monoclonal antibody, PBS Only (Capture) 67528-1-PBS (1H2A6). 100 µg. Concentration 1 mg/mL.

Detection antibody: STK32A Monoclonal antibody, PBS Only (Capture/Detector) 67528-2-PBS (1B8C7). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50267-1, STK32A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67528-1-PBS. Detection antibody: 67528-2-PBS. Standard:Ag10942. Range: 1.563-200 ng/mL.