

For Research Use Only

# ZCCHC3 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50305-2

## Capture Antibody Information

**Catalog Number:**  
68903-1-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG2a  
**Purification Method:**  
Protein A purification

**Clone ID:**  
1G9G3  
**Reactivity:**  
human  
**GenBank:**  
BC069238  
**Immunogen Catalog Number:**  
Ag34163

**Conjugate:**  
Unconjugated  
**Full name:**  
zinc finger, CCHC domain containing 3  
**Gene ID:**  
85364

## Detection Antibody Information

**Catalog Number:**  
68903-3-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
2A5G3  
**Reactivity:**  
human  
**GenBank:**  
BC069238  
**Immunogen Catalog Number:**  
Ag34163

**Conjugate:**  
Unconjugated  
**Full name:**  
zinc finger, CCHC domain containing 3  
**Gene ID:**  
85364

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.781-200 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50305-2 targets ZCCHC3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.  
Capture antibody: ZCCHC3 Monoclonal antibody, PBS Only (Capture) 68903-1-PBS (1G9G3). 100 µg. Concentration 1 mg/ml.  
Detection antibody: ZCCHC3 Monoclonal antibody, PBS Only (Detector) 68903-3-PBS (2A5G3). 100 µg. Concentration 1 mg/ml.  
Alternative ZCCHC3 matched antibody pairs: MP50305-1  
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.  
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.  
Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50305-2, ZCCHC3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68903-1-PBS. Detection antibody: 68903-3-PBS. Standard:Ag34163. Range: 0.781-200 ng/mL.