

TSLP Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50325-1**

Capture Antibody Information

Catalog Number: 68911-1-PBS	Clone ID: 1H1G7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: thymic stromal lymphopoietin
Isotype: IgG1	GenBank: BC016720	Gene ID: 85480
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag4771	

Detection Antibody Information

Catalog Number: 68911-2-PBS	Clone ID: 4G8B11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: thymic stromal lymphopoietin
Isotype: IgG2a	GenBank: BC016720	Gene ID: 85480
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag4771	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50325-1 targets TSLP in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TSLP Monoclonal antibody, PBS Only (Capture) 68911-1-PBS (1H1G7). 100 μ g. Concentration 1 mg/mL.

Detection antibody: TSLP Monoclonal antibody, PBS Only (Detector) 68911-2-PBS (4G8B11). 100 μ g. Concentration 1 mg/mL.

Alternative TSLP matched antibody pairs: MP01287-1, MP01287-2, MP01287-3, MP01287-4, MP50325-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50325-1, TSLP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68911-1-PBS. Detection antibody: 68911-2-PBS. Standard:Ag4771. Range: 0.098-100 ng/mL.