

# c-MYC Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50347-3**

## Capture Antibody Information

**Catalog Number:**  
67447-6-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
4C7  
**Reactivity:**  
human  
**GenBank:**  
BC000141  
**Immunogen Catalog Number:**  
Ag17519

**Conjugate:**  
Unconjugated  
**Full name:**  
v-myc myelocytomatosis viral oncogene homolog (avian)  
**Gene ID:**  
4609

## Detection Antibody Information

**Catalog Number:**  
67447-7-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
4A11B10  
**Reactivity:**  
human  
**GenBank:**  
BC000141  
**Immunogen Catalog Number:**  
Ag17519

**Conjugate:**  
Unconjugated  
**Full name:**  
v-myc myelocytomatosis viral oncogene homolog (avian)  
**Gene ID:**  
4609

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.781-50 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50347-3 targets c-MYC in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MYC Monoclonal antibody, PBS Only (Capture) 67447-6-PBS (4C7). 100 µg. Concentration 1 mg/mL.

Detection antibody: MYC Monoclonal antibody, PBS Only (Detector) 67447-7-PBS (4A11B10). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP50347-3, MYC/c-MYC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67447-6-PBS. Detection antibody: 67447-7-PBS. Standard:Ag17519. Range: 0.781-50 ng/mL.