

# SIGLEC9 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50367-3**

## Capture Antibody Information

**Catalog Number:**  
68937-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG2b  
**Purification Method:**  
Protein A purification

**Clone ID:**  
1A5C4  
**Reactivity:**  
human  
**GenBank:**  
BC035365  
**Immunogen Catalog Number:**  
Ag4201

**Conjugate:**  
Unconjugated  
**Full name:**  
sialic acid binding Ig-like lectin 9  
**Gene ID:**  
27180

## Detection Antibody Information

**Catalog Number:**  
68937-4-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
1B1C11  
**Reactivity:**  
human  
**GenBank:**  
BC035365  
**Immunogen Catalog Number:**  
Ag4201

**Conjugate:**  
Unconjugated  
**Full name:**  
sialic acid binding Ig-like lectin 9  
**Gene ID:**  
27180

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
1.563-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50367-3 targets SIGLEC9 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Siglec-9 Monoclonal antibody, PBS Only (Capture/Detector) 68937-2-PBS (1A5C4). 100 µg. Concentration 1 mg/mL.

Detection antibody: Siglec-9 Monoclonal antibody, PBS Only (Detector) 68937-4-PBS (1B1C11). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP50367-3, SIGLEC9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68937-2-PBS. Detection antibody: 68937-4-PBS. Standard:Ag4201. Range: 1.563-100 ng/mL.