

TRIM21 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50371-1**

Capture Antibody Information

Catalog Number: 67136-2-PBS	Clone ID: 1A5B1	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: tripartite motif-containing 21
Isotype: IgG1	GenBank: BC010861	Gene ID: 6737
Purification Method: Protein G purification	Immunogen Catalog Number: Ag28377	

Detection Antibody Information

Catalog Number: 67136-3-PBS	Clone ID: 1D11D10	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: tripartite motif-containing 21
Isotype: IgG1	GenBank: BC010861	Gene ID: 6737
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28377	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-50 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	--	--

Product Information

MP50371-1 targets TRIM21 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TRIM21 Monoclonal antibody, PBS Only (Capture) 67136-2-PBS (1A5B1). 100 μ g. Concentration 1 mg/mL.

Detection antibody: TRIM21 Monoclonal antibody, PBS Only (Detector) 67136-3-PBS (1D11D10). 100 μ g. Concentration 1 mg/mL.

Alternative TRIM21 matched antibody pairs: MP00505-1, MP50371-2, MP50371-3, MP50371-4

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50371-1, TRIM21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67136-2-PBS. Detection antibody: 67136-3-PBS. Standard:Ag28377. Range: 0.098-50 ng/mL