

# TRIM21 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50371-2**

## Capture Antibody Information

Catalog Number: 67136-1-PBS	Clone ID: 1E8B9	Conjugate: Unconjugated
Host: Mouse	Reactivity: human, mouse, rat	Full name: tripartite motif-containing 21
Isotype: IgG1	GenBank: BC010861	Gene ID: 6737
Purification Method: Protein A purification	Immunogen Catalog Number: Ag28377	

## Detection Antibody Information

Catalog Number: 67136-3-PBS	Clone ID: 1D11D10	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: tripartite motif-containing 21
Isotype: IgG1	GenBank: BC010861	Gene ID: 6737
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28377	

## Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

## Product Information

MP50371-2 targets TRIM21 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TRIM21 Monoclonal antibody, PBS Only (Capture) 67136-1-PBS (1E8B9). 100 µg. Concentration 1 mg/mL.

Detection antibody: TRIM21 Monoclonal antibody, PBS Only (Detector) 67136-3-PBS (1D11D10). 100 µg. Concentration 1 mg/mL.

Alternative TRIM21 matched antibody pairs: MP00505-1, MP50371-1, MP50371-3, MP50371-4

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP50371-2, TRIM21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67136-1-PBS. Detection antibody: 67136-3-PBS. Standard:Ag28377. Range: 0.781-100 ng/mL.