

PRKCQ Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50373-1

Capture Antibody Information

Catalog Number: 68939-1-PBS	Clone ID: 2D1A2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: protein kinase C, theta
Isotype: IgG1	GenBank: BC101465	Gene ID: 5588
Purification Method: Protein G purification	Immunogen Catalog Number: Ag26331	

Detection Antibody Information

Catalog Number: 68939-2-PBS	Clone ID: 1D1G3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: protein kinase C, theta
Isotype: IgG1	GenBank: BC101465	Gene ID: 5588
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag26331	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-25 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP50373-1 targets PRKCQ in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PRKCQ Monoclonal antibody, PBS Only (Capture) 68939-1-PBS (2D1A2). 100 µg. Concentration 1 mg/mL.

Detection antibody: PRKCQ Monoclonal antibody, PBS Only (Detector) 68939-2-PBS (1D1G3). 100 µg. Concentration 1 mg/mL.

Alternative PRKCQ matched antibody pairs: MP50373-2, MP50373-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50373-1, PRKCC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68939-1-PBS. Detection antibody: 68939-2-PBS. Standard:Ag26331. Range: 0.781-25 ng/mL