

# ADAMTS17 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50434-1**

## Capture Antibody Information

<b>Catalog Number:</b> 68973-1-PBS	<b>Clone ID:</b> 2D5F9	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> ADAM metalloproteinase with thrombospondin type 1 motif, 17
<b>Isotype:</b> IgG2a	<b>GenBank:</b> NM_139057.3	<b>Gene ID:</b> 170691
<b>Purification Method:</b> Protein A Magarose purification	<b>Immunogen Catalog Number:</b> Ag32013	

## Detection Antibody Information

<b>Catalog Number:</b> 68973-2-PBS	<b>Clone ID:</b> 1H4C5	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> ADAM metalloproteinase with thrombospondin type 1 motif, 17
<b>Isotype:</b> IgG1	<b>GenBank:</b> NM_139057.3	<b>Gene ID:</b> 170691
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag32013	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.391-50 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
------------------------------------------------------	---------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------

## Product Information

MP50434-1 targets ADAMTS17 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ADAMTS17 Monoclonal antibody, PBS Only (Capture) 68973-1-PBS (2D5F9). 100 µg. Concentration 1 mg/mL.

Detection antibody: ADAMTS17 Monoclonal antibody, PBS Only (Detector) 68973-2-PBS (1H4C5). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP50434-1, ADAMTS17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68973-1-PBS. Detection antibody: 68973-2-PBS. Standard:Ag32013. Range: 0.391-50 ng/mL.