

For Research Use Only

CEP164 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50474-4

Capture Antibody Information

Catalog Number: 67370-7-PBS	Clone ID: 2E8D8	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: centrosomal protein 164kDa
Isotype: IgG1	GenBank: BC000602	Gene ID: 22897
Purification Method: Protein G purification	Immunogen Catalog Number: Ag17699	

Detection Antibody Information

Catalog Number: 67370-3-PBS	Clone ID: 3F1D6	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: centrosomal protein 164kDa
Isotype: IgG1	GenBank: BC000602	Gene ID: 22897
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag17699	

Applications

Tested Applications: Cytometric bead array	Range: 0.391-6.25 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP50474-4 targets CEP164 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CEP164 Monoclonal antibody, PBS Only (Capture) 67370-7-PBS (2E8D8). 100 µg. Concentration 1 mg/ml.

Detection antibody: CEP164 Monoclonal antibody, PBS Only (Detector) 67370-3-PBS (3F1D6). 100 µg. Concentration 1 mg/ml.

Alternative CEP164 matched antibody pairs: MP50474-1, MP50474-2, MP50474-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50474-4, CEP164 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67370-7-PBS. Detection antibody: 67370-3-PBS. Standard: Ag17699. Range: 0.391-6.25 ng/mL