

IL17RA Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50493-1**

Capture Antibody Information

Catalog Number: 60400-1-PBS	Clone ID: 2H4C2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: interleukin 17 receptor A
Isotype: IgG1	GenBank: BC011624	Gene ID: 23765
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag33881	

Detection Antibody Information

Catalog Number: 60400-2-PBS	Clone ID: 2C8G11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: interleukin 17 receptor A
Isotype: IgG3	GenBank: BC011624	Gene ID: 23765
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag33881	

Applications

Tested Applications: Cytometric bead array	Range: 0.195-25 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP50493-1 targets IL17RA in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IL17RA Monoclonal antibody, PBS Only (Capture) 60400-1-PBS (2H4C2). 100 μ g. Concentration 1 mg/mL.

Detection antibody: IL17RA Monoclonal antibody, PBS Only (Detector) 60400-2-PBS (2C8G11). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50493-1, IL17RA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60400-1-PBS. Detection antibody: 60400-2-PBS. Standard:Ag33881. Range: 0.195-25 ng/mL