

For Research Use Only

# LRPPRC Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50550-1

## Capture Antibody Information

Catalog Number:  
67679-2-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
3C8G4

Reactivity:  
human

GenBank:  
BC050311

Immunogen Catalog Number:  
Ag30261

Conjugate:  
Unconjugated

Full name:  
leucine-rich PPR-motif containing

Gene ID:  
10128

## Detection Antibody Information

Catalog Number:  
67679-3-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
1F7B10

Reactivity:  
human

GenBank:  
BC050311

Immunogen Catalog Number:  
Ag30261

Conjugate:  
Unconjugated

Full name:  
leucine-rich PPR-motif containing

Gene ID:  
10128

## Applications

Tested Applications:  
Cytometric bead array

Range:  
3.125-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50550-1 targets LRPPRC in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LRPPRC Monoclonal antibody, PBS Only (Capture/Detector) 67679-2-PBS (3C8G4). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: LRPPRC Monoclonal antibody, PBS Only (Detector) 67679-3-PBS (1F7B10). 100  $\mu$ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50550-1, LRPPRC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67679-2-PBS. Detection antibody: 67679-3-PBS. Standard: Ag30261. Range: 3.125-100 ng/mL.