

CCL24 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50599-3**

Capture Antibody Information

Catalog Number: 60443-3-PBS	Clone ID: 3D11H11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: chemokine (C-C motif) ligand 24
Isotype: IgG1	GenBank: NM_002991	Gene ID: 6369
Purification Method: Protein G purification	Immunogen Catalog Number: Eg0395	

Detection Antibody Information

Catalog Number: 60443-4-PBS	Clone ID: 3E1A1	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: chemokine (C-C motif) ligand 24
Isotype: IgG1	GenBank: NM_002991	Gene ID: 6369
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Eg0395	

Applications

Tested Applications: Cytometric bead array	Range: 0.195-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
------------------------------------------------------	----------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------

Product Information

MP50599-3 targets CCL24 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CCL24/Eotaxin 2 Monoclonal antibody, PBS Only (Capture) 60443-3-PBS (3D11H11). 100 µg. Concentration 1 mg/mL.

Detection antibody: CCL24/Eotaxin 2 Monoclonal antibody, PBS Only (Detector) 60443-4-PBS (3E1A1). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50599-3, CCL24 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60443-3-PBS. Detection antibody: 60443-4-PBS. Standard:Eg0395. Range: 0.195-100 ng/mL.