

# uPAR/CD87 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: MP50621-1

## Capture Antibody Information

**Catalog Number:**  
60456-1-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
3E5C2  
**Reactivity:**  
human  
**GenBank:**  
NM\_002659.4

**Conjugate:**  
Unconjugated  
**Full name:**  
PLAUR  
**Gene ID:**  
5329

## Detection Antibody Information

**Catalog Number:**  
60456-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
3A1B3  
**Reactivity:**  
human  
**GenBank:**  
NM\_002659.4

**Conjugate:**  
Unconjugated  
**Full name:**  
PLAUR  
**Gene ID:**  
5329

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.098-12.5 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50621-1 targets uPAR/CD87 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: uPAR/CD87 Monoclonal antibody, PBS Only (Capture) 60456-1-PBS (3E5C2). 100 µg. Concentration 1 mg/mL.

Detection antibody: uPAR/CD87 Monoclonal antibody, PBS Only (Detector) 60456-2-PBS (3A1B3). 100 µg. Concentration 1 mg/mL.

Alternative uPAR/CD87 matched antibody pairs: MP00129-1, MP00129-2, MP00129-3, MP00129-4, MP00129-5, MP50621-2, MP50621-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP50621-1, uPAR/CD87 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60456-1-PBS. Detection antibody: 60456-2-PBS. Standard:Eg0901. Range: 0.098-12.5 ng/mL.