

HSPD1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50644-3**

Capture Antibody Information

Catalog Number: 66041-2-PBS	Clone ID: 1D6F11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: heat shock 60kDa protein 1 (chaperonin)
Isotype: IgG1	GenBank: BC003030	Gene ID: 3329
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag7848	

Detection Antibody Information

Catalog Number: 66041-4-PBS	Clone ID: 2E8E2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: heat shock 60kDa protein 1 (chaperonin)
Isotype: IgG1	GenBank: BC003030	Gene ID: 3329
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag7848	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-200 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP50644-3 targets HSPD1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HSPD1 Monoclonal antibody, PBS Only (Capture/Detector) 66041-2-PBS (1D6F11). 100 μ g. Concentration 1 mg/mL.

Detection antibody: HSPD1 Monoclonal antibody, PBS Only (Detector) 66041-4-PBS (2E8E2). 100 μ g. Concentration 1 mg/mL.

Alternative HSPD1 matched antibody pairs: MP50644-1, MP50644-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50644-3, HSPD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66041-2-PBS. Detection antibody: 66041-4-PBS. Standard:Ag7848. Range: 1.563-200 ng/mL.