

MTHFD2L Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50702-3**

Capture Antibody Information

Catalog Number: 60504-3-PBS	Clone ID: 1A6B11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: methylentetrahydrofolate dehydrogenase (NADP+ dependent) 2-like
Isotype: IgG1	GenBank: BC032771	Gene ID: 441024
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag11470	

Detection Antibody Information

Catalog Number: 60504-2-PBS	Clone ID: 1C1D7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: methylentetrahydrofolate dehydrogenase (NADP+ dependent) 2-like
Isotype: IgG1	GenBank: BC032771	Gene ID: 441024
Purification Method: Protein G purification	Immunogen Catalog Number: Ag11470	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50702-3 targets MTHFD2L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MTHFD2L Monoclonal antibody, PBS Only (Capture/Detector) 60504-3-PBS (1A6B11). 100 µg. Concentration 1 mg/mL.

Detection antibody: MTHFD2L Monoclonal antibody, PBS Only (Detector) 60504-2-PBS (1C1D7). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50702-3, MTHFD2L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60504-3-PBS. Detection antibody: 60504-2-PBS. Standard:Ag11470. Range: 1.563-100 ng/mL