

# SMCR7L Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50714-4**

## Capture Antibody Information

<b>Catalog Number:</b> 67808-2-PBS	<b>Clone ID:</b> 3A7C9	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> Smith-Magenis syndrome chromosome region, candidate 7-like
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC002587	<b>Gene ID:</b> 54471
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag13775	

## Detection Antibody Information

<b>Catalog Number:</b> 67808-6-PBS	<b>Clone ID:</b> 1C12C3	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> Smith-Magenis syndrome chromosome region, candidate 7-like
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC002587	<b>Gene ID:</b> 54471
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag13775	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.391-100 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

## Product Information

MP50714-4 targets SMCR7L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SMCR7L Monoclonal antibody, PBS Only (Capture) 67808-2-PBS (3A7C9). 100 µg. Concentration 1 mg/mL.

Detection antibody: SMCR7L Monoclonal antibody, PBS Only (Detector) 67808-6-PBS (1C12C3). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP50714-4, SMCRTL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67808-2-PBS. Detection antibody: 67808-6-PBS. Standard:Ag13775. Range: 0.391-100 ng/mL