

AEG-1/MTDH Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50732-1**

Capture Antibody Information

Catalog Number: 68592-2-PBS	Clone ID: 1A8G2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: metadherin
Isotype: IgG1	GenBank: BC045642	Gene ID: 92140
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag5033	

Detection Antibody Information

Catalog Number: 68592-3-PBS	Clone ID: 2C1A2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: metadherin
Isotype: IgG1	GenBank: BC045642	Gene ID: 92140
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag5033	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP50732-1 targets AEG-1/MTDH in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: AEG-1 Monoclonal antibody, PBS Only (Capture) 68592-2-PBS (1A8G2). 100 µg. Concentration 1 mg/mL.

Detection antibody: AEG-1 Monoclonal antibody, PBS Only (Detector) 68592-3-PBS (2C1A2). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50732-1, AEG-1/MTDH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68592-2-PBS. Detection antibody: 68592-3-PBS. Standard:Ag5033. Range: 0.098-100 ng/mL