

# AEG-1/MTDH Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50732-4**

## Capture Antibody Information

Catalog Number: 68592-6-PBS	Clone ID: 1F5E2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: metadherin
Isotype: IgG1	GenBank: BC045642	Gene ID: 92140
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag5033	

## Detection Antibody Information

Catalog Number: 68592-4-PBS	Clone ID: 2B11F8	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: metadherin
Isotype: IgG1	GenBank: BC045642	Gene ID: 92140
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag5033	

## Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

## Product Information

MP50732-4 targets AEG-1/MTDH in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: AEG-1 Monoclonal antibody, PBS Only (Capture) 68592-6-PBS (1F5E2). 100 µg. Concentration 1 mg/mL.

Detection antibody: AEG-1 Monoclonal antibody, PBS Only (Detector) 68592-4-PBS (2B11F8). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage buffer:  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP50732-4, AEG-1/MTDH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68592-6-PBS. Detection antibody: 68592-4-PBS. Standard:Ag5033. Range: 0.781-100 ng/mL