

SSBP1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50750-2**

Capture Antibody Information

Catalog Number:
67582-4-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1A12H1
Reactivity:
human
GenBank:
BC000895
Immunogen Catalog Number:
Ag29976

Conjugate:
Unconjugated
Full name:
single-stranded DNA binding protein 1
Gene ID:
6742

Detection Antibody Information

Catalog Number:
67582-5-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1D3E11
Reactivity:
human
GenBank:
BC000895
Immunogen Catalog Number:
Ag29976

Conjugate:
Unconjugated
Full name:
single-stranded DNA binding protein 1
Gene ID:
6742

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50750-2 targets SSBP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SSBP1 Monoclonal antibody, PBS Only (Capture) 67582-4-PBS (1A12H1). 100 μ g. Concentration 1 mg/mL.

Detection antibody: SSBP1 Monoclonal antibody, PBS Only (Capture/Detector) 67582-5-PBS (1D3E11). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50750-2, SSBP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67582-4-PBS. Detection antibody: 67582-5-PBS. Standard:Ag29976. Range: 0.391-100 ng/mL.