

ASGR2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50757-2**

Capture Antibody Information

Catalog Number: 66831-2-PBS	Clone ID: 1G7F12	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: asialoglycoprotein receptor 2
Isotype: IgG2b	GenBank: BC017251	Gene ID: 433
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag2037	

Detection Antibody Information

Catalog Number: 66831-4-PBS	Clone ID: 2G10C1	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: asialoglycoprotein receptor 2
Isotype: IgG1	GenBank: BC017251	Gene ID: 433
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag2037	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-25 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	--	--

Product Information

MP50757-2 targets ASGR2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ASGR2 Monoclonal antibody, PBS Only (Capture) 66831-2-PBS (1G7F12). 100 µg. Concentration 1 mg/mL.

Detection antibody: ASGR2 Monoclonal antibody, PBS Only (Detector) 66831-4-PBS (2G10C1). 100 µg. Concentration 1 mg/mL.

Alternative ASGR2 matched antibody pairs: MP50757-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50757-2, ASGR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66831-2-PBS. Detection antibody: 66831-4-PBS. Standard:Ag2037. Range: 0.781-25 ng/mL