

LRMDA Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50759-1

Capture Antibody Information

Catalog Number:
60537-1-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A Magarose purification

Clone ID:
1F10D11
Reactivity:
human
GenBank:
NM_032024
Immunogen Catalog Number:
Ag33944

Conjugate:
Unconjugated
Full name:
chromosome 10 open reading frame 11
Gene ID:
83938

Detection Antibody Information

Catalog Number:
60537-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2B12B4
Reactivity:
human
GenBank:
NM_032024
Immunogen Catalog Number:
Ag33944

Conjugate:
Unconjugated
Full name:
chromosome 10 open reading frame 11
Gene ID:
83938

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50759-1 targets LRMDA in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LRMDA Monoclonal antibody, PBS Only (Capture) 60537-1-PBS (1F10D11). 100 μ g. Concentration 1 mg/mL.

Detection antibody: LRMDA Monoclonal antibody, PBS Only (Capture/Detector) 60537-2-PBS (2B12B4). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50759-1, LRMDA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60537-1-PBS. Detection antibody: 60537-2-PBS. Standard:Ag33944. Range: 0.098-100 ng/mL.