

# SLC39A6/ZIP6 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50771-2**

## Capture Antibody Information

<b>Catalog Number:</b> 60546-3-PBS	<b>Clone ID:</b> 1G6E12	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> solute carrier family 39 (zinc transporter), member 6
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC039498	<b>Gene ID:</b> 25800
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag5676	

## Detection Antibody Information

<b>Catalog Number:</b> 60546-4-PBS	<b>Clone ID:</b> 3A3D10	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> solute carrier family 39 (zinc transporter), member 6
<b>Isotype:</b> IgG3	<b>GenBank:</b> BC039498	<b>Gene ID:</b> 25800
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag5676	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.098-12.5 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

## Product Information

MP50771-2 targets SLC39A6/ZIP6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC39A6 Monoclonal antibody, PBS Only (Capture) 60546-3-PBS (1G6E12). 100 µg. Concentration 1 mg/mL.

Detection antibody: SLC39A6 Monoclonal antibody, PBS Only (Detector) 60546-4-PBS (3A3D10). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP50771-2, SLC39A6/ZIP6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60546-3-PBS. Detection antibody: 60546-4-PBS. Standard:Ag5676. Range: 0.098-12.5 ng/mL.