

# STAR Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50772-3**

## Capture Antibody Information

**Catalog Number:**  
67130-3-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
1D5B1  
**Reactivity:**  
human  
**GenBank:**  
BC010550  
**Immunogen Catalog Number:**  
Ag28698

**Conjugate:**  
Unconjugated  
**Full name:**  
steroidogenic acute regulatory protein  
**Gene ID:**  
6770

## Detection Antibody Information

**Catalog Number:**  
67130-5-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
1B9A1  
**Reactivity:**  
human  
**GenBank:**  
BC010550  
**Immunogen Catalog Number:**  
Ag28698

**Conjugate:**  
Unconjugated  
**Full name:**  
steroidogenic acute regulatory protein  
**Gene ID:**  
6770

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
6.25-200 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50772-3 targets STAR in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: STAR Monoclonal antibody, PBS Only (Capture/Detector) 67130-3-PBS (1D5B1). 100 µg. Concentration 1 mg/mL.

Detection antibody: STAR Monoclonal antibody, PBS Only (Detector) 67130-5-PBS (1B9A1). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP50772-3, STAR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67130-3-PBS. Detection antibody: 67130-5-PBS. Standard:Ag28698. Range: 6.25-200 ng/mL