

CYP17A1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: MP50779-1

Capture Antibody Information

Catalog Number: 60554-1-PBS	Clone ID: 1A5F2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: cytochrome P450, family 17, subfamily A, polypeptide 1
Isotype: IgG1	GenBank: BC062997	Gene ID: 1586
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag5998	

Detection Antibody Information

Catalog Number: 60554-2-PBS	Clone ID: 1E4C12	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: cytochrome P450, family 17, subfamily A, polypeptide 1
Isotype: IgG1	GenBank: BC062997	Gene ID: 1586
Purification Method: Protein G purification	Immunogen Catalog Number: Ag5998	

Applications

Tested Applications: Cytometric bead array	Range: 0.195-6.25 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP50779-1 targets CYP17A1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CYP17A1 Monoclonal antibody, PBS Only (Capture) 60554-1-PBS (1A5F2). 100 µg. Concentration 1 mg/mL.

Detection antibody: CYP17A1 Monoclonal antibody, PBS Only (Detector) 60554-2-PBS (1E4C12). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50779-1, CYP17A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60554-1-PBS. Detection antibody: 60554-2-PBS. Standard:Ag5998. Range: 0.195-6.25 ng/mL