

# Integrin alpha-7 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: MP50795-1

## Capture Antibody Information

**Catalog Number:**  
60564-1-PBS

**Host:**  
Mouse

**Isotype:**  
IgG2a

**Purification Method:**  
Protein A Magarose purification

**Clone ID:**  
2D6A2

**Reactivity:**  
human

**GenBank:**  
BC050280

**Immunogen Catalog Number:**  
Ag19777

**Conjugate:**  
Unconjugated

**Full name:**  
integrin, alpha 7

**Gene ID:**  
3679

## Detection Antibody Information

**Catalog Number:**  
60564-2-PBS

**Host:**  
Mouse

**Isotype:**  
IgG1

**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
1D11H5

**Reactivity:**  
human

**GenBank:**  
BC050280

**Immunogen Catalog Number:**  
Ag19777

**Conjugate:**  
Unconjugated

**Full name:**  
integrin, alpha 7

**Gene ID:**  
3679

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.098-6.25 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50795-1 targets Integrin alpha-7 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Integrin alpha-7 Monoclonal antibody, PBS Only (Capture) 60564-1-PBS (2D6A2). 100 µg. Concentration 1 mg/mL.

Detection antibody: Integrin alpha-7 Monoclonal antibody, PBS Only (Detector) 60564-2-PBS (1D11H5). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP50795-1, Integrin alpha-7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60564-1-PBS. Detection antibody: 60564-2-PBS. Standard:Ag19777. Range: 0.098-6.25 ng/mL